

Worry-free pipetting is finally here with the only tip that clips.

Eliminate the frustration of banging tips on your pipette


With a light touch, each tip is clipped securely on the pipette.

Reduce the risk of loose or leaking tips

When Thermo Scientific™ ClipTip™ pipette tips are on, they are sealed and will not loosen. The ClipTip pipetting system remains sealed until you decide to eject tips, even in the most comprehensive applications, including multiple tip touches and mixing cycles.

Enable consistent and reproducible pipetting from user to user

Higher-quality results and more efficient research.



ClipTip Pipette Tips

- Versatile product range with 12.5-1250 µL 96-format tips and 12.5-125µl 384-format tips.
- Several packaging options from racked to environmentally friendly and space saving reload towers.
- All ClipTip tips are certified free of RNase, DNase, DNA, ATP and endotoxin contamination.
- Low retention technology minimizes liquid retention, and enable maximum sample recovery. Available in ClipTip 20-1250 96-format tips.
- Convenient color coding makes it easy to find the compatible F1-ClipTip or E1-ClipTip pipettes.

Thermo Scientific™ F1-ClipTip™ Manual Pipettes

- Revolutionary Interlocking ClipTip technology provides a secure seal on every channel.
- Extremely light tip attachment and ejection forces to help reduce strain injuries.
- Effortless and secure volume adjustment with light and soundless volume adjustment mechanism and volume lock.
- Light-weight handle design with adjustable finger rest.
- Complete range of single, 8- and 12-channel pipettes within volume range of 1-1000 µL.

Thermo Scientific™ E1-ClipTip™ Electronic Pipettes

- Revolutionary Interlocking ClipTip technology provides a secure seal on every channel.
- Adjustable tip-spacing combined with the electronic tip ejection, and index finger operation ensures more ergonomic and efficient pipetting experience.
- An intuitive personalized user-interface makes even the most complex protocols easier, and efficient to complete.
- Versatile product range from single, 8-, 12- and 16-channel to 6-12-channel adjustable tip spacing pipettes.

Save 90% more time with E1-ClipTip Equalizer pipettes!

Save over 90% more time, and 90% less repetitions. Transfer samples between any labware format, and save time with E1-ClipTip Electronic Adjustable Tip Spacing models.

Equalizer model selector guide

96-Format Equalizer Pipettes					
Channels	Max. spacing	Transfer to/from	↔	Transfer from/to	Min. spacing
6-channel	19.8 mm	<ul style="list-style-type: none"> • Test tube racks • 24-well microplates 	↔	<ul style="list-style-type: none"> • 96-well microplates • Deepwell blocks 	9.0 mm
8-channel	14.2 mm	<ul style="list-style-type: none"> • Microcentrifuge tube racks • 48-well microplates 	↔	<ul style="list-style-type: none"> • 96-well microplates • Deepwell blocks 	9.0 mm
384-Format Equalizer Pipettes					
Channels	Max. spacing	Transfer to/from	↔	Transfer from/to	Min. spacing
8-channel	14.2 mm	<ul style="list-style-type: none"> • 96-well microplates • 48-well microplates • Deepwell blocks • Microcentrifuge tube racks 	↔	<ul style="list-style-type: none"> • 384-well microplates • 384 PCR plates • Agarose Gels 	4.5 mm
12-channel	9.0 mm	<ul style="list-style-type: none"> • 96-well microplates • Deepwell blocks 	↔	<ul style="list-style-type: none"> • 384-well microplates • 384 PCR plates • Agarose Gels 	4.5 mm

PROTECTING YOUR INVESTMENT

We provide a complete portfolio of customized, integrated, and innovative services and support solutions designed to help you improve productivity, reduce total cost of ownership, and ensure compliance throughout your laboratory — from instrument and equipment acquisition to disposition.

- Preventive Maintenance
- Calibration
- Onsite Depot

unity lab services
by Thermo Fisher Scientific

thermoscientific

ДИАМ
современная лаборатория



Sealed in security. Feel the difference.

Thermo Scientific ClipTip Pipetting Systems

Selection Guide

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru



